
# Tutorial

# Gemelli

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dic. '21

This work is licensed under a Creative Commons Attribution 4.0 International License.

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
**Main Operation Frequency for 1H** 63.87 MHz

**Frequency range 1H** 63.57 MHz - 64.18 MHz

**Frequency range Multi Nuclei** Option not available

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
Gemelli

The screenshot shows the Philips MRI software interface. The top menu bar includes Patients, Examination, Review, Analysis, System, and Help. The patient information section displays Registration ID, Date of Birth, and Gender (Male). The examination details show 'RM ENCEFALO IN..., 31-10-2016'. The sequence list is as follows:

Sequence	Checkmark	Orientation
1,1 Survey	✓	
2,1 T1W_3D_TFE	✓	SAG(1)
3,1 3D_Brain_FLAIR	✓	*SAG
4,1 T2W_TSE_HR	✓	axi
5,1 DWI	✓	axi
6,1 SWIp_HR	✓	axi
7,1 T2W_TSE_COR	✓	cor
8,1 T2W_TSE_SAG	✓	SAG
9,1 3D_PRESS_144	✓	spet 3D

The '3D\_PRESS\_144' sequence is selected. The right side of the interface shows a sagittal MRI scan of the brain. A context menu is open over the '3D\_PRESS\_144' sequence, with 'SpectroView' highlighted. Other menu options include QFlow, IViewBold, Advanced Properties..., Copy (Ctrl+C), Paste (Ctrl+V), Expand All, and Save Protocol to Text File... A green arrow points to the 'SpectroView' option.

# SPETTROSCOPIA – MULTIVOXEL

The screenshot displays the Philips SpectroView software interface. The top menu bar includes 'Patients', 'Examination', 'Review', 'Analysis', 'System', and 'Help'. The left sidebar shows patient information: 'Registration ID:', 'Date of Birth:', 'Gender: Male', and a dropdown menu for 'RM ENCEFALO IN..., 31-10-2016'. Below this is a table of sequences for '\* Pediatrico HR'.

Sequence ID	Sequence Name	Orientation
1,1	Survey	
2,1	T1W_3D_TFE	SAG(1)
3,1	3D_Brain_FLAIR	<del>SAG</del>
4,1	T2W_TSE_HR	axi
5,1	DWI	axi
6,1	SWIp_HR	axi
7,1	T2W_TSE_COR	cor
8,1	T2W_TSE_SAG	SAG
9,1	3D_PRESS_144	spet 3D

The right pane shows 'SpectroView - 3D\_PRESS\_144'. An 'Unknown Anatomy' dialog box is open, displaying the message: 'Anatomy is not yet defined for this series. Please specify the anatomy for this study.' The dialog includes a dropdown menu with 'Brain' selected and an 'OK' button.

# SPETTROSCOPIA – MULTIVOXEL

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
...

# SPETTROSCOPIA – MULTIVOXEL

Selezionare:  
«Select  
Script»


# SPETTROSCOPIA – MULTIVOXEL

Please Select a Script

Name	Description	Nucleus	Anatomy	Field	TE
<unnamed>-20161008-110842	Based on Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	long
<unnamed>-20161011-121114	Based on Long-TE Brain (1H CSI or TSI, modu...	1H	Brain	1.5T	long
<unnamed>-20161013-113913	Based on Long-TE Brain (1H CSI or TSI, modu...	1H	Brain	1.5T	long
<unnamed>-20161014-091615	Based on Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.5T	long
<unnamed>-20161014-091708	Based on Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.5T	long
<unnamed>-20161014-091739	Based on Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.5T	long
<unnamed>-20161031-112103	Based on Long-TE Brain (1H CSI or TSI, modu...	1H	Brain	1.5T	long
BasicProc_csi	Basic Processing (1H CSI)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	short
<b>LongTeBrain_csi</b>	<b>Long-TE Brain (1H CSI or TSI, modulus spectra)</b>	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	long
LongTeBrain_csi-real	Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	long

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Name:

Description:

Anatomy:  Supported Field Strengths:

Nucleus:   1.0T  1.5T  3.0T  7.0T

Echo Time:  Short  Long

Processing Step:

Spectrum Phase Adjustment  
 Initial Baseline Subtraction  
 Shift Peak Frequency  
 Select Peaks  
 Peak Fitting  
 Correct for DSA filter  
 Generate Maps  
 Results Table  
 Graph Display

Ac  
Ala  
Cho  
Cr  
Cr2  
Glx  
H2O  
Lac  
lip  
ml  
NAA  
sl

Note: assuming 1.5T peaks!

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# SPETTROSCOPIA – MULTIVOXEL

**Script Parameters**

Name: <unnamed> -20161031-113526

Description: Based on Long-TE Brain (1H CSI or TSI, modulus spectra)

Anatomy: Brain Supported Field Strengths:

Nucleus: 1H  1.0T  1.5T  3.0T  7.0T

Echo Time:  Short  Long

Processing Step: Select Peaks

Spectrum Phase Adjustment  
 Initial Baseline Subtraction  
 Shift Peak Frequency  
 Select Peaks  
 Peak Fitting  
 Correct for DSA filter  
 Generate Maps  
 Results Table  
 Graph Display

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