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
MRI – PROPELLER - MULTIVANE

Gemelli

This work is licensed under a Creative Commons Attribution 4.0 International License.

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
@Live Expert session

Propeller Imaging

Maria Cammarata
MR Clinical Education Specialist
IMI

14 Marzo 2014

ore 17.30

Durata : 1 ora

e-mail:

maria.cammarata@ge.com

PROPELLER IMAGING

Propeller DWI **Propeller DW Scan 1**

PROPELLER coil-HEAD, mode-1, table-0.0

APPLICATOR PARAMETERS

Phase Position:	Superior	Phase:	180.0 deg
Poloid Slice:	Head First	Frequency:	120
TE (sec):	102.0	Slice Thickness (mm):	5.0
TR (sec):	5000.0	Slice-spacing (mm):	3.0
PTV (cm):	23.0	Cartis Frequency:	
Headwidth (cm):	31.250	Auto Stack Length:	20
TE (sec) (2000 Hz):	100.0	Auto Stack Interval (sec):	40.0

APPLICATION OPTIONS

Full Data Collection

Auto

Graphical

Acquire: _____

Speed: _____

DIFFUSION PROFILES

All

Auto


Command

L1
 R1
 AP

L1
 R1
 AP

Diffusion Coefficient:
 Diffusion Time:
 Save Raw Data

Save Series Scan Time: _____ WAcq: D Max Shims/Acq: 20 Total Shims: 0 TR: 617
 -dddt Est. SAR: _____ Avg. Coi SAR: _____ Peak SAR: _____

Periodically Rotated Overlapping Parallel Lines with Enhanced Reconstruction

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



PROPELLER imaging significantly reduces the effects of motion artifact in routine T1, T2, T2FLAIR, T1 FLAIR, and PD scanning. DW PROPELLER significantly improves image quality in the vicinity of bone/tissue or air/tissue interfaces, or around tissue/metal interfaces prone to creating susceptibility artifacts.

PROPELLER K-Space filling:

In comparison to FSE, PROPELLER acquires k-space distinctly different. PROPELLER fills k-space in a unique way. Rather than going at it line by line, k-space is filled with an arrangement of "blades." The blades' radial trajectory removes structured motion artifact, and redundant sampling enables the reduction of bulk patient motion artifacts.

Proprietary Training Material Property of GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1

gen. '24

21

Protocol Optimization Guidelines for PROPELLER:

- **No Phase Wrap:** Care should be taken to select the proper **No Phase Wrap factor** to reduce streaking or wrap artifacts. With the introduction of AIR Recon in DV28 software, NPW factors can be reduced for time savings. As a general guideline, No Phase Wrap factors are provided below:
 - o Brain and Spine 1.5 NPW
 - o MSK 1.5 NPW (for Transmit Receive coils) or 2 NPW (for Receive only coils)
 - o Pelvis 2 to 2.5 NPW

Note! To gain SNR in your PROPELLER with AIR Recon DL sequences, it is recommended to INCREASE NEX rather than increase No Phase Wrap factors to avoid potential for streak artifacts

Proprietary Training Material Property of GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1

gen. '24

22

Protocol Optimization Guidelines for PROPELLER:

Bandwidth:

Bandwidth selection is especially important when referring to PROPELLER protocols utilizing AIR Recon DL. Selecting an appropriate receive bandwidth (rBW) minimizes chemical shift artifacts and reduces echo spacing. Follow the indications below:

- ✓ 3T: up to 125 kHz with AIR Recon DL
- ✓ 1.5T: up to 83 kHz with AIR Recon DL

PROPELLER IMAGING

Protocol Optimization Guidelines for PROPELLER: Acceleration:

With AIR Recon DL, you can apply acceleration techniques for a time savings.

- Minimum ARC of 2 provides a time savings without sacrificing image quality

Pictured Below: A) Coronal T2 PROPELLER knee sequence using rBW of 62 kHz, ARC 3, 2:12 demonstrating reduced chemical shift artifact (circled in green) and a time savings; B) Same sequence using 31 kHz, ARC 3, 2:59 demonstrating increased chemical shift artifact (circled in red)

Protocol Optimization Guidelines for PROPELLER:

Echo Train Length (ETL):

For T2 FLAIR and STIR PROPELLER sequences, the ETL selection is grayed out to reduce streak artifacts, especially when using AIR Recon DL. To effectively balance your ETL, adjust your TE using the equation below:

- $(TE/\text{Echo Spacing}) \times 2$

TE (Echo Time): In general, the TE times for PROPELLER sequences are higher than FSE sequences:

- For PD contrast: +10 ms (~20%)
- For T2 contrast: +20 ms (~20%)

Proprietary Training Material Property of GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1

gen. '24

25

Protocol Optimization Guidelines for PROPELLER:

Note!

For **1.5T systems***: If you have upgraded your system from MR29.1 or lower to **MR30** with AIR Recon DL, you may see your TE values blank out. This is expected behavior due to the removal of User CV5 Editable TE. If your values have blanked, you can select an appropriate value using the guidelines below. o Note!

For **3T systems***: If you have upgraded your system from MR29 or lower to MR30 with AIR Recon DL for PROPELLER, you may see your TE values blank out. This is expected behavior due to the removal of User CV5 Editable TE. If your values have blanked, you can select an appropriate value using the guidelines below.

Proprietary Training Material Property of GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1

gen. '24

26

Protocol Optimization Guidelines for PROPELLER:

TE to 100 ms for T2 with ETL greater than or equal to 16

TE to 45ms for PD with ETL between 8 and 16

TE to 20ms for T1 with ETL equal to or less than 8

TE to 45ms for T1 FLAIR with ETL between 8 and 16

TE to 100 ms for T2 FLAIR with ETL greater than or equal to 16

** For systems on MR30.1 or higher, your protocol should automatically populate with an appropriate value*

Protocol Optimization Guidelines for PROPELLER:

Harmonize: Harmonize corrects for gross motion in PROPELLER sequences, and is recommended for HEAD ONLY imaging •

ChemSat Duo: Recommended and useful for MSK off center imaging in 1.5T Systems

NEX: PROPELLER with AIRTM Recon DL is locked to a MINIMUM 1.5 NEX

Note!

To gain SNR in your PROPELLER with AIR Recon DL sequences, it is recommended to INCREASE NEX rather than increase No Phase Wrap factors to avoid potential for streak artifacts.

PROPELLER IMAGING

Protocol Optimization Guidelines for PROPELLER:

PROPELLER T1 Protocol Review:

Neuro: T1 optimization ON, NEX \geq 2, ETL $>$ 5 Flip Angle 111

MSK: T1 optimization OFF, ETL \geq 5, Flip Angle 111

Spine: T1 optimization ON, NEX \geq 2, ETL \geq 5

General Considerations:

- The GE Protocol Library has optimized AIR™ Recon DL protocols for all anatomies. Use these as a starting guideline.
- Minimum 1.5 NEX
- Increase rBW

Pictured: *SIGNA Premier* **Coronal T1 PROPELLER**
Knee 384x384; NPW 1.5, NEX 1.5; 83 kHz; 2min
13sec

PROPELLER IMAGING

Protocol Optimization Guidelines for PROPELLER:


PROPELLER T1 FLAIR Protocol Review:

Spine: Flip Angle 130

General Considerations:

- The GE Protocol Library has optimized AIR™ Recon DL protocols for all anatomies. Use these as a starting guideline.
- Minimum 1.5 NEX
- Increase rBW

Pictured: *SIGNA Architect Sagittal T1 FLAIR PROPELLER Cervical Spine*; 268x268; NEX 2; 83.3 kHz; 3min 32sec

Proprietary Training Material Property of GE HealthCare © 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1

PROPELLER IMAGING


Protocol Optimization Guidelines for PROPELLER:

PROPELLER T2 Protocol Review:

Spine: ETL roughly TE/Echo Spacing *2

General Considerations:

- The GE Protocol Library has optimized AIR™ Recon DL protocols for all anatomies. Use these as a starting guideline.
- Minimum 1.5 NEX
- Increase rBW
- TE should be roughly 20% higher than your FSE protocols to achieve the same contrast

Pictured: (Left) *SIGNA Premier Axial T2 PROPELLER Brain* 512x512; NEX 1.5; 91 kHz; 2min 52sec (Right) *SIGNA Premier Axial T2 PROPELLER Orbits* 384x384; NPW 1.3; NEX 1.5; 83 kHz; 1 min 24 sec

PROPELLER IMAGING

Protocol Optimization Guidelines for PROPELLER:

PROPELLER T2 FLAIR Protocol Review:

General Considerations:

- The GE Protocol Library has optimized AIR™ Recon DL protocols for all anatomies. Use these as a starting guideline.
- Minimum 1.5 NEX
- Increase rBW

Pictured: *SIGNA Architect Axial T2 FLAIR PROPELLER* Brain 300x300; NEX 1.76; 83kHz; 4min 41sec

Proprietary Training Material Property of GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1

gen. '24

32

PROPELLER IMAGING

Protocol Optimization Guidelines for PROPELLER:


PROPELLER STIR Protocol Review:

For STIR PROPELLER and FSE, new optimal TI may have changed with STIR enhancements in MR30 or higher software for **3T Systems only**. Historically, TI for PROPELLER STIR sequences was 180, and is now 210 ms.

General Considerations:

- The GE Protocol Library has optimized AIR™ Recon DL protocols for all anatomies. Use these as a starting guideline.
- Minimum 1.5 NEX
- Increase rBW

Pictured: *SIGNA Premier Sagittal STIR PROPELLER*
Cspine 288x288; NPW 1.5, NEX 2; 62.5kHz; 2min
40sec

Proprietary Training Material Property of GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1

gen. '24

33

PROPELLER IMAGING

Protocol Optimization Guidelines for PROPELLER:

PROPELLER PD Protocol Review:


MSK: TE should be roughly 10% higher than your FSE protocols to achieve the same contrast

General Considerations:

- The GE Protocol Library has optimized AIR™ Recon DL protocols for all anatomies. Use these as a starting guideline.
- Minimum 1.5 NEX
- Increase rBW

Pictured: *SIGNA Premier Axial PD*

PROPELLER Wrist 540x540; NPW 1.4, NEX 1.5; 91 kHz; 4min 16sec

Proprietary Training Material Property of GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1

PROPELLER IMAGING

Protocol Optimization Guidelines for PROPELLER:



DW PROPELLER Protocol Review:

Neuro: Duo (For Silent DW PROPELLER, use dwprep)

General Considerations:

- The GE Protocol Library has optimized AIR™ Recon DL protocols for all anatomies. Use these as a starting guideline.
- Minimum 1.5 NEX
- Increase rBW

Pictured: (Top) *SIGNA Architect Axial DWI PROPELLER Brain 128x128; NPW 1.1; 83.3 kHz; 1 min 52 sec;* (Bottom) *SIGNA Architect Cor DWI PROPELLER Brain 128x128; NPW 1.2; 83.3 kHz; 2 min 7 sec*

Proprietary Training Material Property of GE HealthCare

© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.

DOC 2815354 Rev 1


gen. '24

35

PROPELLER IMAGING

Artifacts:

Streaking artifact is a known limitation of PROPELLER. Increasing NEX from 1.5 to 3 significantly reduces this artifact.


Pictured above: (Left) Axial T1 PROPELLER Brain using 3 NEX with no streaking artifact (circled in green); (Right) Axial T1 PROPELLER Brain using 1.5 NEX with apparent streaking (circled in red)

Proprietary Training Material Property of GE HealthCare © 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license. DOC 2815354 Rev 1

Protocol Optimization Guidelines for PROPELLER:

Zipper artifacts present along the phase encoding direction for each blade, similar to FSE. Zipper artifacts in each blade are rotated according to the corresponding blade angle and they appear similar to streaking artifacts in the final image. Zipper artifacts are generally more intense than streaking due to motion. To minimize zipper artifacts, ensure that there are no RF leaks in the room due to an open door, faulty light bulb, etc.

Pictured below: Axial T2 PROPELLER Pelvis zipper artifacts are present in each blade along the phase encoding direction

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MultiVane XD

	Technical characteristics	Additional information
Sequence	Multislice TSE or FFE	Routine clinical contrasts: T2w, FLAIR, all orientations
Main applications	Brain, Abdomen, Spine, HeadNeck	Also available for other anatomies
Field strength	1.5T and 3.0T	
Parallel imaging	Compatible with dS SENSE	Leverages the efficient, high SNR, dS SENSE parallel imaging technology to further reduce scan time to routine levels