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
# AIR RECON DL

## Protocol Optimization For AIR™ Recon DL Recon DL

## Protocol Optimization For AIR™ Recon DL

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## Protocol Optimization For AIR™ Recon DL

### 3T Protocols:

- +/- **42.5 to 62 kHz** ((or higher)  
Increasing bandwidth allows  
for a decreased echo spacing
- BW/pixel value at **350-450**

### 1.5T Protocols:

- +/- **32kHz to +/-42kHz**, is  
sufficient to sharpen images  
without adversely impacting  
remaining SNR.
- BW/pixel value at **300- 350**

## Protocol Optimization For AIR™ Recon DL

Why do I still see motion on my images? AIR™ Recon DL removes noise ringing artifacts from your images but note that it does not remove additional artifacts from your images. In fact, removing the noise from your image may increase the likelihood of noticing other artifacts, such as motion. How can you avoid motion on your images?

- Make sure your patient is comfortable on the scan table, providing pillows and pads as necessary
- Immobilize the anatomy as much as possible
- Use PROPELLER sequences Especially in Cervical Spine Imaging: • Remind your patient to minimize swallowing during scans

## Protocol Optimization For AIR™ Recon DL

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DOC2723319 Rev.01

# AIR RECON DL

## Level up your image quality with AIR™ Recon DL

AIR™ Recon DL is a pioneering, deep-learning based reconstruction algorithm that improves SNR and image sharpness, enabling shorter scan times. It improves image quality at the foundational level by making use of the raw data to remove image noise and ringing.

Increases productivity by enabling shorter scan times

Removes image noise and ringing by leveraging raw data

Delivers sharper and clearer TrueFidelity™ MR images


Enables you to set your preferred SNR improvement level



## AIR™ Recon DL MR image reconstruction

Deep-learning network trained to reconstruct images that are:

- sharper from Intelligent True Resolution
- higher in SNR with a user-selectable improvement level

Raw data




Final image

## Air Recon DL: Cervical Cancer

FSE T2 1.5 NEX **2:47**

FSE T2 1 NEX **:58**  
**Increased Acc**

FSE T2 DL Recon **:58**



↑ Same raw data ↑

## AIR RECON DL – HOW TO USE IT?

PROTOCOL OPTIMIZATION !

## EASE-OF-USE

- AIR Recon DL shows impressive performance in reducing image noise and Gibbs Ringing and improves SNR/CNR
- Allows to push protocols or applications to levels incompatible with conventional reconstruction – in any anatomy.
- Avoids the tradeoff between resolution, scan time and image quality. Gives you the „free lunch“.
- Most significant development in MR technology over the last several years.

	Conventional	AIR Recon DL	% Decrease
Shoulder	12:49	<b>6:38</b>	50%
LSP	11:04	<b>6:44</b>	40%
CSP	14:06	<b>8:24</b>	40%
TSP	9:38	<b>6:27</b>	35%
Knee	13:27	<b>8:35</b>	35%
Prostate	17:39	<b>12:50</b>	30%
Hip	17:58	<b>9:23</b>	45%
Wrist	13:47	<b>7:50</b>	45%
Hand	13:07	<b>5:47</b>	50%
Ankle	12:22	<b>7:50</b>	40%
Foot	11:54	<b>7:27</b>	40%
Female Pelvis	17:54	<b>12:09</b>	35%
Abdomen	13:02	<b>7:01</b>	40%

## AIR Recon DL

- Elevates image quality
  - Reduces image noise
  - Recognizes and mitigates artifacts
  - Encourages use of higher matrices (resolution)
  - Enables techniques that reduce blur & inherently sharpens
  - Reduces motion artifact
- Substantially improves patient experience and throughput from 33-50%

## Protocol Optimization For AIR™ Recon DL

ZIP1024

ZIP = zero-filling  
that enhances  
apparent image  
resolution

**Ax DWI ALL b1000**    GRx    0:44

Details    Diffusion    Acceleration    Advanced

Scan Plane: **Oblique**    Freq. Dir: **R/L**

Freq. FOV: **24.0**    Auto TR: **8816.0**

Phase FOV: **1.00**    # Slices: **40**

Slice Thickness: **3.00**

Spacing: **0.3**

	R/L	A/P	S/I	Max # Slices	40
Start	<b>L2.7</b>	<b>A17.0</b>	<b>S55.1</b>	# of Acqs	<b>1</b>
End	<b>R8.1</b>	<b>A10.6</b>	<b>S80.7</b>	Rel. SNR(S)	<b>70</b>
Act End	<b>R7.5</b>	<b>A11.0</b>	<b>S73.1</b>	Acq Voxel Size	<b>1.5x1.5x3.0</b>
Chem SAT: <b>Fat</b>				Echo Spacing	<b>0.9</b>

Intensity Filter: **None**

Save Original:

3D Geometry Correction:

Real Time Field Adjustment:

Real Time Center Frequency:

Minimum TE: 77.4

Maximum TE: 169

WB-SAR: 0.17    Head-SAR: 0.52    B<sub>1</sub>-RMS: 0.97µT    Mode: First    dB/dt: First

**Ax DWI ALL b1000**    GRx    0:44    Details    **Diffusion**    Acceleration    Advanced

Scan Plane: Oblique    Freq. Dir: R/L    # b-values: 1

Freq. FOV: 24.0    Auto TR: 8816.0    NEX for T2: 1.00

Phase FOV: 1.00    # Slices: 40    Diffusion Direction: ALL

Slice Thickness: 3.00    # of Diffusion Directions: 3

Spacing: 0.3    # of T2 Images: 1

R/L    A/P    S/I    Max # Slices: 40    # Synthetic b-values: 0

Start L2.7    A17.0    S55.1    # of Acqs: 1

End R8.1    A10.6    S80.7    Rel. SNR: 70

Act End R7.5, A11.0, S73.1    Acq Voxel Size: 1.5x1.5x3.0

Chem SAT: Fat    Echo Spacing: 0.9

Recon All Images:     Optimize TE:     Dual Spin Echo:

#	b-value	NEX
1	1000.0	1.0

## AIR™ Recon DL – Protocol Basics

The Test

You need to get 6 correct answers to pass this test.

Once you view the replay and pass the test you will be able to download and Save or Print the CE Certificate.

1. ARDL is a deep-learning based reconstruction algorithm that improves SNR and image sharpness, enabling shorter scan times.

Select a Choice

-

-

TRUE

2. FALSE

2. Which statement is/are correct regarding ARDL?

Select a Choice

-

-

It is a T2 weighted sequence

3. Can only be used in AXIAL images

Enables you to set your preferred SNR improvement level

All of the above

3. ARDL removes image noise and ringing artifacts.

Select a Choice

-

-

TRUE

4. FALSE

## AIR™ Recon DL – Protocol Basics

4. ARDL is compatible with all anatomies and various pulse sequences.

Select a Choice

-

TRUE

FALSE

5. Which parameters are more commonly adjusted to optimize your protocols for ARDL?

Select a Choice

-

Matrix, BW, NEX, Recon DL strength

6. TR, TE, Frequency direction

Scan Plane, number of slices

None of the above

6. ZIP1024 is recommended when using ARDL

Select a Choice

-

TRUE

FALSE

7. Which parameter should be adjusted on DWI if using ARDL?

Select a Choice

-

NEX & MATRIX

8. NEX

MATRIX

SPACING

## AIR™ Recon DL – Protocol Basics

8. What is the recommended Bandwidth at 1.5T and 3T when using ARDL?

Select a Choice

-

+/- 10 to 16kHz for 1.5T and +/- 16 to 31.5kHz for 3T

+/- 32 to 42kHz for 1.5T and +/- 42.5 to 62kHz for 3T

None of the above

All of the above

**TiP**<sup>®</sup> (Training in Partnership)

GE Healthcare certifies that

[Redacted]

has satisfactorily completed

AIR Recon DL - Protocol Basics On Demand

and has been awarded 1

Continuing Education credits as designated by

ASRT Department of Education – Category A credits

July 21, 2022

08/01/2025

Margie Travis

GE Representative

WIZ00320139

Date Reference Number

Activity/Course Expiration License Number

Please retain this certificate for your own personal record of continuing education. Certificate may be replaced only by GE Healthcare TiP Applications if damaged or lost.

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This work is licensed under a Creative Commons Attribution 4.0 International License.