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



```
NEC MultiSync P330W
sdev on /dev type tmpfs (rw)
loop: module loaded
devpts on /dev/pts type devpts (rw,mde=0520,gid=5)
kjournald starting. Commit interval 5 seconds
EXT3 FS on sda3, internal journal
EXT3-fs: mounted filesystem with ordered data mode.
/dev/sda3 on /export/home type ext3 (rw,acl,user_xattr)
kjournald starting. Commit interval 5 seconds
EXT3 FS on sda5, internal journal
EXT3-fs: mounted filesystem with ordered data mode.
/dev/sda5 on /export/home/signa/research type ext3 (rw,acl,user_xattr)
kjournald starting. Commit interval 5 seconds
EXT3 FS on sda6, internal journal
EXT3-fs: mounted filesystem with ordered data mode.
/dev/sda6 on /export/home/snap type ext3 (rw,acl,user_xattr)
kjournald starting. Commit interval 5 seconds
EXT3 FS on md0, internal journal
EXT3-fs: mounted filesystem with ordered data mode.
/dev/md0 on /export/home1 type ext3 (rw,smntime)
Retry device configuration:
Setting current systemctl status from /etc/systemd.conf
net.ipv4.icmp_echo_ignore_broadcasts = 0
net.ipv4.conf.all.rp_filter = 1
fs.notify_max_user_watches = 65536
net.ipv4.conf.default.promote_secondaries = 1
net.ipv4.conf.all.promote_secondaries = 1
kernel.core_pattern = /export/home/service/log/core.%s.%t.jp

Enabling syn flood protection
Disabling IP forwarding

Activating remaining swap-devices in /etc/fstab...
fuse init (API version 7.9)
Loading fuse module
Mounting fuse control filesystem
Creating /var/log/boot.msg
Mounting securityfs on /sys/kernel/security
Loading AppArmor profiles

Setting up hostname 'ix-nr'
Setting up loopback interface lo
lo IP address: 127.0.0.1/8
lo IP address: 127.0.0.2/8

Starting ldconfig
System Boot Control: The system has been
Skipped features:
System Boot Control: Running /etc/init.d/boot.local
INIT: Entering runlevel: 5
Boot logging started on /dev/tty1(/dev/console) at Tue Jan 19 07:01:59 2016
Master Resource Control: previous runlevel: N, switching to runlevel:
acpid: starting up

Starting acpid
Initializing random number generator
acpid: 2 rules loaded

Starting srodUSB
Checking/Updating CPU microcode
Starting syslog services
Starting D-Bus daemon
Starting timeZoneStartupScript
Loading CPUFreq modules
Starting INL daemon_
```

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


Inserire login  
e password...

**SIGNA™ Premier**

© 2020 General Electric Company - All Rights Reserved

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




NEC MultiSync P241W

Patient Record: Name: Find: Search/Refresh Last Refreshed Time: 19 Jan, 07:28:43

Worklist T... Scheduled... Name & Patient ID Accession SPS Description Modal... Referring P... Scheduled... Status Local Protocol Refs... Caution

**Patient**

Last Name First Name  
Name

Patient ID:

lbs Kgs  
Weight

DD-MM-YYYY  
Date Of Birth:

Year Month Weeks Days  
Age:

Sex:

DD-MM-YYYY  
Scheduled Date:

\* Required

**Exam**

Accession:

Exam Description:

SPS Description:

Procedure ID:

Last Name First Name  
Operator

Last Name First Name  
Radiologist

Last Name First Name  
Referring Physician

Protocol: Enter Number

Favorite Protocols: Select Favorite

Show All Protocols...

**Other Information**

This is patient sensitive information, please adhere to patient confidentiality rules. Please verify the information using either the patients medical record or directly with the patient to guarantee its accuracy.

Allergies: NONE

Pre-Med: NONE

Pregnancy Status: No Entry

History:

SAR Display

	1.0 Sec	5 Min	Room °C
Avg W/kg			
Limit W/kg	6.0	10.0	
Act W/kg	0.0	0.0	

7:10 19 January 21%

THE TABLE IS NOT DOCKED.

Scanner  
not assigned to a patient

**Patient Record**  
 In Progress    Details  
 Completed

**Find:**    Last:    First:

**Search/Refresh**  
 Last Refreshed Time: 21 Dec, 07:16:31

SPS Description	Schedul...	Worklist...	Name ▲	Patient ID	Accession	Status	Modality	Schedu...	Local Pr...	Referr...	Recom...	Caution
[Empty Table Area]												

**Patient**  
 Last Name:     First Name:   
 Patient ID:   
 Weight:  lb     kg  
 Date Of Birth:  DD-MM-YYYY  
 Age:  year(s)    month(s)    week(s)    day(s)  
 Sex:

**Exam**  
 Accession:   
 Exam Description:   
 SPS Description:   
 Procedure ID:   
 Operator:  Last Name    First Name  
 Radiologist:  Last Name    First Name  
 Referring Physician:  Last Name    First Name  
 Protocol:  Enter Number  
 Favorite Protocols:  Select Favorite

**Other Information**  
 This is patient sensitive information, please adhere to patient confidentiality rules. Please verify the information using either the patients medical record or directly with the patient to guarantee its accuracy.  
 Other Patient ID:   
 Current Location:   
 Special Needs:   
 Requested Contrast Agent:   
 Allergies: NONE  
 Pre-Med: NONE  
 Pregnancy Status:  No Entry  
 History:

Full View     Report Cursor     Auto W/L  

## GE Healthcare

**SAR Display**

Avg W/kg	10 Sec	6 Min	Room °C
Limit W/kg:	6.0	10.0	18.0
Act W/kg:	0.0	0.0	18.0

7:18  
21 December  
47%

THE TABLE IS NOT AT SCANNING HEIGHT.

The image displays five screenshots of an MRI console interface, numbered 1 through 5:

- 1:** Patient information screen. Fields include Patient Name (Doe Jane), Patient ID (1234567), DOB (11-12-2000), Exam, Patient Weight (123 lbs.), Patient Position (diagram), and Landmark (diagram). A timer shows 00:00.
- 2:** Vital signs monitoring screen. Shows ECG, SpO2 (80%), and Respiration (15). A timer shows 01:47.
- 3:** Scan parameters screen. Shows a large '12 mm' and a timer of 01:47.
- 4:** Coil selection screen. Lists 'Port 1 Flex Large 16ch' and 'Port E Spine Posterior 60ch'. Includes a 'Notes' section with a 'Port Compatibility Key' (16 Ch. or less, 16 to 32 Ch.) and a T/R indicator. A timer shows 01:47.
- 5:** Fan/Light Control screen. Shows controls for fan and light with a timer of 01:47.

## MR Systems Ingenia

Ingenia (2017-09-19)

Philips Medical Systems Nederland B.V.  
Veenpluis 4-6, 5684 PC Best, The Netherlands.

Logon

Name

Password


OK

Cancel

PHILIPS

Healthcare

PHILIPS

Healthcare

Patients Examination Review Analysis System Help

Registration ID: Plan  
Date of Birth: Review  
Gender:

00:00:00

Add new scan item...

PHILIPS  
Healthcare

Voxel Rel SNR TE TR Accept SAR PNS Cancel

Summary

FOV x mm x mm  
Voxel x mm x mm  
Matrix x x slices  
Gap  Default mm


NSA  
Fat saturation  SPIR  
Low SAR mode

Info Assistance AutoView

Start Scan Stop Scan

2023  
BEST HOSPITALS  
Newsweek  
1° OSPEDALE IN ITALIA  
statista  
POLICLINICO UNIVERSITARIO A. GEMELLI

PHILIPS

Healthcare

### New Examination

**Patient**

Patient name: \*

Registration ID: \*

Date of birth: \*  dd-mm-yyyy

Age:

Gender: \*

Patient weight: \*  kg

**Examination**

Exam name:

Accession number:

Examination date:   21-12-2020

Referring Physician:

Performing Physician:

Study Comments:

Allowed SAR mode:

**Patient conditions**

Pregnant \*

Implant \*

Medical alerts:

Allergies:

Start Date	Start TI...	Patient Name	Registration ID	Accession N...	Date Of Bl...	Gender	Procedure Step
------------	-------------	--------------	-----------------	----------------	---------------	--------	----------------

No exams scheduled.

fenix  
the healthcare solution by elco

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Refreshed Time : 19 Jan , 07:09:43

Search Data

Fetch the RIS and the locally scheduled worklists.

Local Protocol

Reco...

Caution

GE Healthcare

Patient Record

View:  In Progress  Details  Completed

Name: [Dropdown] Last: [Text] First: [Text] Search Data [Button]

Worklist T... Scheduled ... Name ▲ Patient ID Accession SPS Description Modal... Referring P... Scheduled ... Status Local Protocol Reco... Caution

08:00:00 [Text] [Text] [Text] RM cranio con MDC in anestesia MR 19-01-2016

### Patient

Last Name: FONTANA First Name: GIOELE

\* Patient ID: 28973899

Lbs: [Text] Kgs: [Text]

\* Weight: [Text] [Text]

DD-MM-YYYY

Date Of Birth: 11-04-2013

Year: [Text] Month: [Text] Weeks: [Text] Days: [Text]

Age: 2

Sex: M

DD-MM-YYYY

Scheduled Date: 19-01-2016

\* Required

### Exam

Accession: IU116GG3DXCX

Exam Description: RM cranio con MI

SPS Description: RM cranio con MI

Procedure ID: REQ2016001752

Operator: Last Name: [Text] First Name: [Text]

Radiologist: Last Name: [Text] First Name: [Text]

Referring Physician: Last Name: [Text] First Name: [Text]

Protocol: Enter Number

Favorite Protocols: Select Favorite

Show All Protocols...

### Other Information

This is patient sensitive information, please adhere to patient confidentiality rules. Please verify the information using either the patient's medical record or directly with the patient to guarantee its accuracy.

Allergies: NONE

Pre-Med: NONE



Pregnancy Status: No Entry

History:

[Text Area]

Start Exam [Button]


GE Healthcare

<b>Radiologist</b>	<b>Last Name</b>	<b>First Name</b>
	<input type="text"/>	<input type="text"/>
<b>Referring Physician</b>	<b>Last Name</b>	<b>First Name</b>
	<input type="text"/>	<input type="text"/>
<b>Protocol:</b>	<input type="text" value="Enter Number"/>	
<b>Favorite Protocols:</b>	<input type="text" value="Select Favorite"/> ▼	
	<input type="button" value="Show All Protocols..."/>	


**History:**

GE Healthcare

GE Healthcare

The screenshot displays a software interface for configuring an MRI protocol. At the top, there is a 'Service' button. Below it, a 'Filter List by:' section contains several checkboxes: 'Standard', 'Express', 'Advance', 'Vascular', 'Contrast', and 'Other'. A central window titled 'Sorted by Alphabetical ▲' lists various MRI sequences, including 'LOCALIZER', 'ASSET Cal', 'Sag T2 FSE', 'Ax T2 FSE', 'Ax T1 FSE', 'Ax T2 FLAIR 4 min', 'Cor T2 FSE', 'Ax DWI TETRA B1000', 'AX 3D SWAN', 'AX 3D FSPGR + C', 'Cor T1 FSE + c fs', 'Sag T1 FSE + c fs', 'Ax T1 FSE + c fs', 'Perfusion', 'DTI (25 Directions)', 'Ax T2 FSE spettro', 'PROBE-SV 144', and 'PROBE-SI 144 Multi voxel'. A 'Multi Protocol Basket' window is empty on the right. At the bottom, a 'Protocol Description' field shows '11.15-\*\*\*Enc Neoplasie Pediatrico 2015' and 'Protocol Description: Protocollo Ottimizzato'. Below that, 'Tasks' are listed as 'LOCALIZER' and 'ASSET Cal'. At the very bottom, patient information is visible: 'ID: 28973899', 'Age: 2Y', 'Lbs/Kgs', and 'Description: RM cranio con MDC in anestesia'. 'Accept' and 'Cancel' buttons are at the bottom right.

**Exam dB/dt and SAR Limits**

Select the desired Operating Mode:

**dB/dt:**

- Normal dB/dt
- First Level dB/dt


Entering First Level for dB/dt.  
Operator/Physician confirms awareness of potential risks and accepts the responsibility.

**SAR:**

- Normal SAR
- First Level SAR

Entering First Level for SAR.  
Operator/Physician confirms awareness of potential risks and accepts the responsibility.

Follow proper patient padding and positioning instructions to prevent patient warming.

**Accept** **Normal Mode**

GE Healthcare

GE Healthcare

5016 \*\*\*Enc Neoplasie Pediatrico 2... 72:05

Exam Edit End


Task	Series Data
InRx	LOCALIZER 00:32
	ASSET Cal 00:12
	Sag T2 FSE 03:25
	Ax T2 FSE 03:40
	Ax T1 FSE 02:33
	Ax T2 FLAIR 4 mm 04:01
	Cor T2 FSE 04:11
	Ax DWI TETRA B1000 03:21
	AX 3D SWAN 04:06
	AX 3D FSPGR +C 04:19
	Cor T1 FSE +c fs 04:02
	Sag T1 FSE +c fs 03:01
	Ax T1 FSE +c fs 03:33

Setup Add Task ▶ Run

View

Auto Start  Auto Scan Scan

Save Rx

Coil...  BHRBRAIN Head First, Supine

Imaging Options... 2D Spin Echo Seq, EDR, Fast, SS

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
### General

### Push Nodes

Patient Position	Supine
Patient Orientation	Head First
Laterality:	Unpaired
Anatomy:	Head
Anatomic Region:	Brain
Table Usage	Use
Heart Rate (beats/min)	60
Align Overlap (mm)	30
GeoLink Propagation	No
Geometries	AXX
Disengage posterior coil	No

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When prompted, click on the reason for the system shutdown: ***Daily, Service, Other*** or ***Cancel*** to exit the shutdown of your system.

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This work is licensed under a Creative Commons Attribution 4.0 International License.