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# TC – PROCEDURE

# Gemelli

This work is licensed under a Creative Commons Attribution 4.0 International License.

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
## Powering on the CT system

### Power Distribution Unit (PDU)


The **PDU** supplies power to various parts of the system, including components for the stand, table, and operator console. The **PDU** supplies **power** to various parts of the system, including components for the stand, table, and operator console.

On the front of the PDU there are **controls** that are used to **indicate** whether power is connected to the unit, a **button to turn on / off** the power of the gantry and table (marked in red) and an **emergency stop button**.

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
The **Reset power to the Gantry button** may be flashing if a reset is required. (After full power shutdown or if emergency buttons were pressed)  
Remember to press the **Reset button on the gantry control if it is flashing.**  
See the **attention messages** that may appear in the **Tube heating section.**

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Est. max Z location CTDIvol: 19.94 mGy  
 Projected series DLP: 789.46 mGy-cm  
 Accumulated exam DLP: 0.00 mGy-cm

Series Description:

Add Group    Split Current Group    Delete Selected Group    More Info    Smart Prep Rx    ODM

Gating 65 BPM    S. Assist    Prior    Next

Images	Scan Type	Start Location	End Location	No. of Images	Thick Speed	Interval (mm)	Tilt	SFOV	kV	mA	Total Exposure Time	Prep Group (s)	ISD (s)	Breath Hold (s)	Breathe Time (s)	Voice Lights Timer	Clive Duration (s)
1-191	Helical Full 0.35 s	388.750	1148.750	191	1.25 20.62 0.516:1	1.250	30.0	Large Body	120	50	4.73	3.4	1.7	N	N	2 T	0.3

Cancel    Prior Series    Next Series    Create New Series    Repeat Series    Delete Series    Priority Recon    Image Check    Auto Scan    Accept

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### Smart Prep

**On** **Dynamic Transition** **Show Localizer**

**AutoVoice Pre-Message**

Monitor Location	mA	Monitoring Delay	Monitoring ISD	Enhancement Threshold	Diagnostic Delay	Auto Minimum Delay
SD.00	40	10.0	3.0	50	3.2	

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AutoVoice Pre-Message 2

Auto Minimum Delay

Monitor Location	mA	Monitoring Delay	Monitoring ISD	Enhancement Threshold	Diagnostic Delay
50.00	40	10.0	2.5	50	3.1

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Monitor Location	mA	Monitoring Delay	Monitoring ISD	Enhancement Threshold	Diagnostic Delay
50.00	40	10.0	2.5	50	5.4

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**Name:** ID:178    **Protocol:** 5.1 Chest Routine    **Exam:** 104    **Series:** 1

**Anatomical Reference:** SN

**Patient Orientation:** Feet First

**Patient Position:** Supine

**Filming:** AutoFilm Setup, Camera Laser Camera

**Auto Store**    **Auto Transfer PACS**    **Dose Report Auto Transfer**    **Dose SR Report Auto Transfer**

**Copy Pt.Orient. Pt.Position Anat.Ref.**

**Gantry Display Patient Info:** On

**Series Description:** Chest native 0.625mm

**Real Time Information**

**Patient Handling**

**Scanning**

**Injecting**

**Delay Timer**    **Monitor Scan Time** 0.00

**Scan Progress**

Seconds: 0 10 20 30 40 50 60

**Baseline Phase**

**Monitor Phase including 10 sec Monitoring Delay**

**Scan Phase**

Images	Type	kV	mA	Time	Thk	Tilt	Start	End
1-1	Axial	120	40	0.40	5.0	\$0.0	\$0.000	\$0.000
2-41	Axial	120	40	16.00	5.0	\$0.0	\$0.000	\$0.000
1-241	Helical	120	100	6.62	0.625	\$0.0	\$0.000	I300.000

**End Exam**    **Next Series**    **Monitor Phase**    **Scan Phase**    **Priority Recon**    **Image Check**    **Repeat Last Group**

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

AutoFilrn Setup

Camera  
None

Auto Store    Auto Transfer    Dose Report Auto Transfer    Dose SR Report Auto Transfer

Xtream Injector    Show Localizer

**Dose Information**    Setup

Attention: Unrecognized tube - Dose not validated by GE

Images	CTDIvol mGy (NV)	DLP mGy-cm	Dose Eff. %	Phantom cm
1-241	29.96 (N)	997.82	92.60	Head 16

Est. max Z location CTDIvol: 29.96 mGy  
Projected series DLP: 997.82 mGy-cm  
Accumulated exam DLP: 0.00 mGy-cm

Series Description

Add Group	Split Current Group	Delete Selected Group	More Info	Smart Prep Rx	ODM	Gating	S. Assist	Prior	Next								
Images	Scan Type	Start Location	End Location	No. of Images	Thick Speed	Interval (mm)	Gantry Tilt	SFOV	kV	mA	Total Exposure Time	Prep Group (s)	ISD (s)	Breath Hold (s)	Breathe Time (s)	Voice Lights Timer	Cine Duration (s)
1-241	Helical Full 0.5 s	30.000	1300.000	241	0.625 39.37 0.984:1	1.250	30.0	Head	120	500 40.00 <sup>~</sup> ODH	4.40	0.0	1.3	N	N	N	2.0


Cancel
Prior Series
Next Series
Create New Series
Repeat Series
Delete Series
Priority Recon
Image Check
Auto Scan
Accept

New Patient
 Emergency Patient
 Patient Schedule
 Protocol Management
 Retro Recon
 Recon Mgmt
 Daily Prep
 Scanner Utilities

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**Standard** (C) Cardiac compatible

### Soft (C)

**Soft#** Results in finer texture and similar noise compared to **Soft**. Only for head SFOV.

**Stnd (C)** For Routine Soft tissue imaging

**Stnd#** Results in finer texture and similar noise compared to **Stnd**. Only for head SFOV

**Detail (C) (G)** where hybrid tissue detail and bone edges are important, Detail is to standard as standard to soft

**Lung (C)** for interstitial lung pathology.

### Bone (C)

### Bone Plus

**Ultra:** For improved inner ear detail

**Edge** for small bone work in the head, as well as high resolution scans

### Edge Plus\*

**Chest** for mediastinum and lung detail studies. It provides soft tissue resolution and contrast when viewing the images in a soft tissue/mediastinal W/L, and high resolution of the lung tissue when viewing the images in a lung W/L.

Optimize not Needed

Gating


Select the desired Reconstruction Algorithm.

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