
Tutorial

Gemelli  1964
2024

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mag. '24

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This work is licensed under a Creative Commons Attribution 4.0 International License.

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

Gemelli

Spettroscopia Multi-Voxel

The screenshot shows the Philips SpectroView software interface. On the left, there is a patient information panel with fields for Registration ID, Date of Birth, and Gender (Male). Below this is a list of MRI sequences:

Sequence ID	Sequence Name	Orientation
1,1	Survey	
2,1	T1W_3D_TFE	SAG(1)
3,1	3D_Brain_FLAIR	SAG
4,1	T2W_TSE_HR	axi
5,1	DWI	axi
6,1	SWIp_HR	axi
7,1	T2W_TSE_COR	cor
8,1	T2W_TSE_SAG	SAG
9,1	3D_PRESS_144	spet 3D

An 'Unknown Anatomy' dialog box is overlaid on the right side of the screen. The dialog contains the following text:

Unknown Anatomy
Anatomy is not yet defined for this series.

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

Spettroscopia Multi-Voxel

Please Select a Script

Name	Description	Nucleus	Anatomy	Field	TE
<unnamed>-20161008-110842	Based on Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	long
<unnamed>-20161011-121114	Based on Long-TE Brain (1H CSI or TSI, modu...	1H	Brain	1.5T	long
<unnamed>-20161013-113913	Based on Long-TE Brain (1H CSI or TSI, modu...	1H	Brain	1.5T	long
<unnamed>-20161014-091615	Based on Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.5T	long
<unnamed>-20161014-091708	Based on Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.5T	long
<unnamed>-20161014-091739	Based on Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.5T	long
<unnamed>-20161031-112103	Based on Long-TE Brain (1H CSI or TSI, modu...	1H	Brain	1.5T	long
BasicProc_csi	Basic Processing (1H CSI)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	short
LongTeBrain_csi	Long-TE Brain (1H CSI or TSI, modulus spectra)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	long
LongTeBrain_csi-real	Long-TE Brain (1H CSI, real spectra)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	long

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Name:

Description:

Anatomy: Supported Field Strengths:

Nucleus: 1.0T 1.5T 3.0T 7.0T

Echo Time: Short Long

Processing Step:

- Spectrum Phase Adjustment
- Initial Baseline Subtraction
- Shift Peak Frequency
- Select Peaks
- Peak Fitting
- Correct for DSA filter
- Generate Maps
- Results Table
- Graph Display

All None

Ac

Ala

Cho

Cr

Cr2

Glx

H2O

Lac

lip

ml

NAA

sl

Run OK Cancel

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...

Spettroscopia Multi-Voxel

Script Parameters

Name:

Description:

Anatomy: Supported Field Strengths:

Nucleus: 1.0T 1.5T 3.0T 7.0T

Echo Time: Short Long

Processing Step:

Spectrum Phase Adjustment
 Initial Baseline Subtraction
 Shift Peak Frequency
 Select Peaks
 Peak Fitting
 Correct for DSA filter
 Generate Maps
 Results Table
 Graph Display

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