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Protocol Optimization For AIR™ Recon DL

3T Protocols:

- **+/- 42.5 to 62 kHz** ((or higher)
Increasing bandwidth allows
for a decreased echo spacing
- BW/pixel value at **350-450**

1.5T Protocols:

- **+/- 32kHz to +/-42kHz**, is
sufficient to sharpen images
without adversely impacting
remaining SNR.
- BW/pixel value at **300- 350**

Protocol Optimization For AIR™ Recon DL

Why do I still see motion on my images? AIR™ Recon DL removes noise ringing artifacts from your images but note that it does not remove additional artifacts from your images. In fact, removing the noise from your image may increase the likelihood of noticing other artifacts, such as motion. How can you avoid motion on your images?

- Make sure your patient is comfortable on the scan table, providing pillows and pads as necessary
- Immobilize the anatomy as much as possible
- Use PROPELLER sequences Especially in Cervical Spine Imaging: • Remind your patient to minimize swallowing during scans

Protocol Optimization For AIR™ Recon DL

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[gehealthcare.com/mr](https://www.gehealthcare.com/mr)

DOC2723319 Rev.01

Level up your image quality with AIR™ Recon DL

AIR™ Recon DL is a pioneering, deep-learning based reconstruction algorithm that improves SNR and image sharpness, enabling shorter scan times. It improves image quality at the foundational level by making use of the raw data to remove image noise and ringing.

Increases productivity by enabling shorter scan times

Removes image noise and ringing by leveraging raw data

Delivers sharper and clearer TrueFidelity™ MR images


Enables you to set your preferred SNR improvement level



AIR™ Recon DL MR image reconstruction

Deep-learning network trained to reconstruct images that are:

- sharper from Intelligent True Resolution
- higher in SNR with a user-selectable improvement level


Raw data


Final image

Air Recon DL: Cervical Cancer


FSE T2 1.5 NEX **2:47**

FSE T2 1 NEX **:58**
Increased Acc

FSE T2 DL Recon **:58**



Same raw data

AIR RECON DL – HOW TO USE IT?

PROTOCOL OPTIMIZATION !

EASE-OF-USE


- AIR Recon DL shows impressive performance in reducing image noise and Gibbs Ringing and improves SNR/CNR
- Allows to push protocols or applications to levels incompatible with conventional reconstruction – in any anatomy.
- Avoids the tradeoff between resolution, scan time and image quality. Gives you the „free lunch“.
- Most significant development in MR technology over the last several years.

	Conventional	AIR Recon DL	% Decrease
Shoulder	12:49	6:38	50%
LSP	11:04	6:44	40%
CSP	14:06	8:24	40%
TSP	9:38	6:27	35%
Knee	13:27	8:35	35%
Prostate	17:39	12:50	30%
Hip	17:58	9:23	45%
Wrist	13:47	7:50	45%
Hand	13:07	5:47	50%
Ankle	12:22	7:50	40%
Foot	11:54	7:27	40%
Female Pelvis	17:54	12:09	35%
Abdomen	13:02	7:01	40%

AIR Recon DL

- Elevates image quality
 - Reduces image noise
 - Recognizes and mitigates artifacts
 - Encourages use of higher matrices (resolution)
 - Enables techniques that reduce blur & inherently sharpens
 - Reduces motion artifact
- Substantially improves patient experience and throughput from 33-50%

Protocol Optimization For AIR™ Recon DL

ZIP = zero-filling
that enhances
apparent image
resolution

AIR RECON DL

Ax DWI ALL b1000 **GRx** **0:44** **Details** **Diffusion** **Acceleration** **Advanced**

Scan Plane: **Oblique** Freq. Dir: **R/L** # of TE(s) per Scan: **1.0** Frequency: **160**

Freq. FOV: **24.0** Auto TR: **8816.0** Num Shots: **1** Phase: **160**

Phase FOV: **1.00** # Slices: **40** TE: **minimum** Bandwidth: **258.0**

Slice Thickness: **3.00** Recon DL Strength: **High** Excitation Mode: **Select**

Spacing: **0.3** Intensity Correction: **PURE** Shim: **Auto**

R/L A/P S/I Max # Slices: **40** Calibration In Prescan: **On** RF Drive Mode: **Preset**

Start **L2.7** **A17.0** **S55.1** # of Acqs: **1** Intensity Filter: **None** Anatomical Reference: **ACPC**

End **R8.1** **A10.6** **S80.7** Rel. SNR(s): **70** Save Original: ☐ Anatomical Plane: **Axial**

Act End **R7.5, A11.0, S73.1** Acq Voxel Size: **1.5x1.5x3.0** 3D Geometry Correction: ☐ Table Delta: **0.00**

Chem SAT: **Fat** Echo Spacing: **0.9** Real Time Field Adjustment: ☒ Real Time Center Frequency: ☐

☐ Contrast

WB-SAR: 0.17 Head-SAR: 0.52 B₁-RMS: 0.97μT Mode: First dB/dt: First Minimum TE: 77.4 Maximum TE: 169

AIR RECON DL

Ax DWI ALL b1000 **GRx** **0:44** **Details** **Diffusion** **Acceleration** **Advanced**

Scan Plane: **Oblique** Freq. Dir: **R/L** # b-values: **1**

Freq. FOV: **24.0** Auto TR: **8816.0** NEX for T2: **1.00**

Phase FOV: **1.00** # Slices: **40** Diffusion Direction: **ALL**

Slice Thickness: **3.00** # of Diffusion Directions: **3**

Spacing: **0.3** # of T2 Images: **1**

Synthetic b-values: **0**

#	b-value	NEX
1	1000.0	1.0

R/L A/P S/I

Start **L2.7** **A17.0** **I55.1**

End **R8.1** **A10.6** **S80.7**

Act End **R7.5, A11.0, S73.1**

Chem SAT: **Fat**

Max # Slices: **40**

of Acqs: **1**

ReL SNR: **70**

Acq Voxel Size: **1.5x1.5x3.0**

Echo Spacing: **0.9**

☐ Contrast ☐ Recon All Images ☒ Optimize TE ☐ Dual Spin Echo

AIR™ Recon DL – Protocol Basics

The Test

You need to get 6 correct answers to pass this test.

Once you view the replay and pass the test you will be able to download and Save or Print the CE Certificate.

1. ARDL is a deep-learning based reconstruction algorithm that improves SNR and image sharpness, enabling shorter scan times.

Select a Choice

-
- TRUE
- FALSE

2. Which statement is/are correct regarding ARDL?

2. Which statement is/are correct regarding ARDL?

Select a Choice

-
- It is a T2 weighted sequence
- Can only be used in AXIAL images
- Enables you to set your preferred SNR improvement level
- All of the above

3. ARDL removes image noise and ringing artifacts.

Select a Choice

-
- TRUE
- FALSE

AIR™ Recon DL – Protocol Basics

4. ARDL is compatible with all anatomies and various pulse sequences.

Select a Choice

6. ZIP1024 is recommended when using ARDL

Select a Choice

5. Which parameters are more commonly adjusted to optimize your protocols for ARDL?

Select a Choice

6. TR, TE, Frequency direction

7. Which parameter should be adjusted on DWI if using ARDL?

Select a Choice

8. NEX

AIR™ Recon DL – Protocol Basics

8. What is the recommended Bandwidth at 1.5T and 3T when using ARDL?

Select a Choice

-

-
- +/- 10 to 16kHz for 1.5T and +/- 16 to 31.5kHz for 3T
- +/- 32 to 42kHz for 1.5T and +/- 42.5 to 62kHz for 3T
- None of the above
- All of the above

TiP®

(Training in Partnership)

GE Healthcare certifies that

has satisfactorily completed

AIR Recon DL - Protocol Basics On Demand

and has been awarded 1

Continuing Education credits as designated by

ASRT Department of Education – Category A credits

July 21, 2022

WIZ00320139

Date

Reference Number

08/01/2025

Activity/Course Expiration

License Number


Margie Travis

GE Representative

Please retain this certificate for your own personal record of continuing education. Certificate may be replaced only by GE Healthcare TiP Applications if damaged or lost.

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The screenshot displays the 'Advanced' tab of the AIR Recon DL interface. The parameters are organized into two columns:

- Left Column:**
 - # of TE(s) per Scan: 2
 - TE: 1.2
 - TE2: 2.4
 - Flip Angle: 65
 - Recon DL Strength: High** (highlighted with a yellow box)
 - Intensity Correction: PURE
 - Calibration In Prescan: On
 - Intensity Filter: None
 - Save Original: ☐
 - 3D Geometry Correction: ☐
- Right Column:**
 - Frequency: 288
 - Phase: 212
 - Bandwidth: 166.7
 - Excitation Mode: Selecti...
 - Shim: Auto
 - RF Drive Mode: Preset
 - Phase Correct: Off

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This work is licensed under a Creative Commons Attribution 4.0 International License.