
Tutorial

Gemelli

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set. '26

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This work is licensed under a Creative Commons Attribution 4.0 International License.

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set. '26

Unknown Anatomy

Anatomy is not yet defined for this series.

Please specify the anatomy for this study.

Brain


OK

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/view **SpectroView** Print

Current Script: BasicProc_sv

Select Script

Scan Nr. 3, 2 Slice 04/50
T1TFE /M
TE 3.5ms - TR 7.6ms

Please Select a Script

Name	Description	Nucleus	Anatomy	Field	TE
BasicProc_sv	Basic Processing (1H SV)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	short
LongTeBrain_sv	Long-TE Brain (1H SV)	1H	Brain	1.0T, 1.5T, 3.0T, 7.0T	long
SPETTRO	Based on Basic Processing (1...	1H	Brain	1.5T	long

Filter:

- ☒ Nucleus
- ☒ Anatomy
- ☒ Field
- ☒ TE

Script Parameters

Name: SPETTRO

Description: Based on Basic Processing (1H SV)

Anatomy: Brain

Nucleus: 1H

Echo Time: ☐ Short ☐ Long

Supported Field Strengths: ☐ 1.0T ☒ 1.5T ☐ 3.0T ☐ 7.0T

Processing Step: Dual Volume Decode

☒ Dual Volume Decode

☒ Residual Water Subtraction

☐ Time Domain Signal Shifting

☒ Apodization Filtering

☒ Zero-Fill and FT

☐ Spectrum Phase Adjustment


☐ Shift Peak Frequency


☒ Graph Display

There are no parameters for this processing step.

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


16:02 Optimal water suppression final value: 342.8 (deg)
16:01 Performing automatic water suppression optimization...
16:01 FWHM of peak 5.9 Hz
16:00 Iterative shimming on VOI procedure finished
15:59 Performing iterative shimming...

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

```
16:02 Optimal water suppression final value: 342.8 (deg)
16:01 Performing automatic water suppression optimization...
16:01 FWHM of peak 5.9 Hz
16:00 Iterative shimming on VOI procedure finished
15:59 Performing iterative shimming...
```

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




The screenshot shows the Philips SpectroView software interface. On the left, there is a list of MRI sequences under the heading 'Cranio tumori'. The sequence 'SV_PRESS_144' is highlighted. On the right, a dialog box titled 'Unknown Anatomy' is displayed. The dialog box contains the text 'Anatomy is not yet defined for this series.' and 'Please specify the anatomy for this study.' Below this text is a dropdown menu with 'Brain' selected. At the bottom of the dialog box, the 'OK' button is highlighted with a red circle, and a green arrow points to it from below.

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



Script Parameters

Name: <unnamed> -20210126-170148

Description: Based on Long-TE Brain (1H SV)

Anatomy: Brain

Nucleus: 1H

Echo Time: ☐ Short ☒ Long

Processing Step: Initial Baseline Subtraction

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